

PREVENT DRY EYES IN THE AM OF AM BLINK




WITH MODERN LIFE DEMANDING MORE AND MORE TIME SPENT ON OUR PHONES AND COMPUTERS, PEOPLE ARE INCREASINGLY FEELING THE EFFECTS OF DRY EYES. PRESENTING OUR NEW PRODUCT, DRY EYE WAND, WHICH WILL HELP UNBLOCK THE MEIBOMIAN GLANDS. IT LOOKS GOOD ON YOUR SHELVES, IS EASY TO USE.

THE BENEFIT OF USING THIS AGAINST THE TRADITIONAL HEATING COMPRESS METHODS, LIKE TOWELS/NAPKINS IS THAT THE TEMPERATURE WON'T DROP! IT RANGES FROM 37°C TO 45°C. IT IS ALSO SMOOTHER AND LESS HARSH ON YOUR SKIN BECAUSE OF THE DESIGN, AND WON'T CAUSE IRRITATION TO THE SKIN.


SIMPLY APPLY TEAR DROPS, EYE CREAM OR COCONUT OIL ON THE HEAD. PRESS THE EYE WAND SIDEWAYS TO COVER MAX SURFACE AREA, AND SLOWLY GLIDE IT ACROSS FROM INNER TO THE OUTER CORNERS OF YOUR EYELID AND UNDEREYE. LET THE HEAT PENETRATE INTO THE SKIN FOR MAXIMUM BENEFIT. THE VIBRATION ALSO HELPS UNCLOG THE GLANDS FASTER.



Contact us at:
+91 9560087282
sales@exegimed.com



1. Take an appropriate amount of artificial tears or any eye lubricant and spread it evenly on the skin.



2. Long press switch button to turn on device and press switch button to switch the mode, and press "+"/"-" to change temperature.




3. Press the eye wand sideways and gently glide it across your eyelids and undereyes. Let heat penetrate for max benefit.



4. After using each time, long pressing the switch button to turn off, please wipe the massage head with isopropyl alcohol to clean.

it has 3 modes:
vibration+hot compress+red light
hot compress+red light
vibration+hot compress



scan for link to intstructional video